

jetPRIME[®] transfection reagent

Short protocol – siRNA Transfection

DAY 0: Cell seeding

→ Seed cells in **V** mL of serum containing medium according to the table below

Quantities per well, dish or flask

Culture vessel	Number of cells	V = volume of growth medium for cell seeding
24-well	25 000 - 40 000	0.5 mL
12-well	50 000 - 80 000	1 mL
6-well / 35 mm	100 000 - 150 000	2 mL
100 mm / flask 75 cm ²	0.5 x 10 ⁶ - 1 x 10 ⁶	10 mL


DAY 1: Transfection

→ Perform transfection **in the presence of serum**

→ Use **jetPRIME[®] buffer only**

→ Transfect cells at **50% confluency**

Dilute **X** pmoles of siRNA
in **W** µL of jetPRIME[®] buffer
Vortex 10 s and spin down



Quantities per well, dish or flask

Culture vessel	W = volume of jetPRIME [®] buffer	X = amount of siRNA added (10nM)	X = amount of siRNA added (50nM)	Y = volume of jetPRIME [®] reagent
24-well	50 µL	5.5 pmoles (76 ng)	27.5 pmoles (381 ng)	2 µl
12-well	100 µL	11 pmoles (152 ng)	55 pmoles (762 ng)	3 µl
6-well / 35 mm	200 µL	22 pmoles (306 ng)	110 pmoles (1524 ng)	4 µl
100 mm / flask 75 cm ²	500 µL	105 pmoles (1460 ng)	525 pmoles (7274 ng)	20 µL

DAY 2-3: Measure gene expression

See back page for optimization tips

Download complete protocol on <http://www.polyplus-transfection.com/resources/product-literature/>

jetPRIME[®] transfection reagent

Short protocol – Optimization Tips (siRNA)

+ Protocol Optimization

- + Check our online Cell Transfection Database for cell specific protocols at:
<http://www.polyplus-transfection.com/resources/cell-transfection-database/>
- + Test different siRNA concentration ranging from 10 to 50 nM (final concentration).
- + Use cells at 50% confluency at time of transfection.



+ Tips to increase cell viability of sensitive cells

- + Replace medium 4 h after transfection.
- + Check that silencing the target gene does not affect cell viability.

+ Use appropriate controls

- + Positive control: housekeeping gene (GAPDH or HPRT) or fluorescently labeled siRNA.
- + Negative control: mismatch, scramble or non-targeting sequence.

+ Good siRNA Transfection Practices

- + Store appropriately jetPRIME[®] ($5 \pm 3^{\circ}\text{C}$).
- + Ensure that cells have been passaged more than twice and less than 20 times prior to transfection.
- + Discard overconfluent cells.
- + Check the half-lives of the protein and mRNA of interest and measure gene silencing accordingly 24 to 96 h after transfection.
- + Regularly check for mycoplasma contaminations.
- + Serum quality may drastically affect transfection efficiency. When purchasing a new batch of serum or trypsin, check cell viability as well as transfection efficiency.

Note: jetPRIME[®] is also recommended for DNA transfection, virus production and DNA/siRNA cotransfection, please refer to the complete protocol available online at:

<http://www.polyplus-transfection.com/resources/product-literature/>