

	Bead material	Silica			Zirconium					Steel	Garnet with Zr Satellite	Silica and Zirconium
	Bead diameter	0.1mm	0.5mm	1.0mm	0.1mm	0.5mm	1.0mm	1.5mm	3.0mm	2.8mm	0.5/6mm	Mixed
Animal	Adipose							X				
	Artery							X	X			
	Bone									X		
	Brain							X	X			
	Insect - fly								X	X		
	Insect - Tick								X	X		
	Insect - Mosquito						X	X	X			
	Gonad							X	X			
	Hair								X	X	X	
	Heart							X	X			
	Kidney							X				
	Liver											
	Lung							X	X			
	Muscle								X	X		
	Pancreas							X	X			
	Skin								X		X	
Spleen							X	X				
Teeth								X	X			
Microbial	Algae		X			X						
	Bacteria	X	X		X	X						
	Fungi		X	X		X	X					
	Mold		X			X						
	Spores		X	X		X	X					
	Yeast		X	X		X	X					
Plant and Soil	Leaves								X	X		
	Nuts									X		
	Plant, soft tissue								X	X		
	Plant, tough tissue								X	X	X	
	Seeds								X	X		
	Stems								X	X		
	Soil							X				X
	Coral									X		
Fecal	Stool							X				X
	Scat								X	X		X